(EN) Intended use

Eurotrol AlbuTrol® is an assayed albumin control intended for professional use in the verification of the precision and accuracy of the HemoCue® Albumin 201 system. Eurotrol AlbuTrol® is intended for the quantitative determination of albumin.

Summary and Principle

Eurotrol AlbuTrol® is a quality control material available at two levels, each with a known albumin concentration. Eurotrol AlbuTrol® may be used in the periodic verification of the precision and accuracy of the HemoCue® Albumin 201 system. Please refer to local guidelines for recommended frequency of use.

Reagents

Eurotrol AlbuTrol® is provided at two physiologically relevant levels. Each vial contains 1.0 mL of human urine containing albumin.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- Caution: Eurotrol AlbuTrol® is prepared from human urine and contains human source material. Components of the control, which are derived from human source material, have been tested and found to be non-reactive for HBsAg, anti-HCV and anti-HIV. However, no test method can offer complete assurance that products derived from human source material are free from infectious agents.
- This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.

Storage and stability

Eurotrol AlbuTrol® must be stored in the refrigerator at 2–8 °C (36–46 °F). Stored unopened at this temperature, stability is guaranteed until the expiration date as indicated on the vial and the outer box. After opening, Eurotrol AlbuTrol® is stable for 30 days stored at 2–8 °C (36–46 °F).

Procedure

- For optimal performance the product should be treated as a patient sample and instructions for use of the HemoCue® Albumin 201 system should be followed.
- Equilibrate the vial for 15 minutes to a temperature of 15 to 30°C (59–86°F).
- Directly before use, shake the vial by hand for at least 30 seconds to completely mix the contents.
- Hold the vial in an upright position and tap the cap to remove liquid from the dropper tip (Figure 1).
- Open the cap carefully without squeezing the vial (Figure 2).
- Discard the first drop of the quality control liquid.
- Apply a drop on a hydrophobic surface and use it immediately for measurement as indicated in the analyzer user manual.
- After sampling, remove any liquid of the quality control from the colored cap with a clean tissue and recap properly.
- Return to indicated storage conditions after use. The entire procedure has to be repeated every time upon removing the quality control from storage during the entire open-vial time period.

Limitations

- The quality control material should not be used past open-vial time period or expiration date.
- Eurotrol AlbuTrol® is not to be used as a calibrator / standard.
- The ranges of the product are intended only as a guideline and laboratories should determine the ranges based on their own test system and tolerance limits.

Expected values

The mean assay value provided on the enclosed assay sheet is derived from replicate analyses on factory calibrated HemoCue® Albumin 201 systems and are batch specific.

For further information and/or technical assistance, please contact:

	Eurotrol B.V. Keplerlaan 20, 6716 BS Ede, The Netherlands Phone +31 318 695 777 E-mail office@eurotrol.com www.eurotrol.com	Eurotrol Inc. 850 N Black Branch Rd, Elizabethtown, KY USA Phone +1 502 501 1180 E-mail officeUSA@eurotrol.com
---	--	--

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state where the user and/or the patient is established.

Revision history

Revision	Effective Date	Change to previous revision
01	2010-05-18	First revision
02	2011-10-06	Revision of section 'Expected results', change of layout, correction of symbols
03	2017-07-24	Addition of Portuguese translation, editing of symbols, new document ID
04	2018-02-28	Revised contact information, editing of symbol description, new document ID
05	2019-01-23	Addition of Spanish translation, complete revision of IFU including new layout

All trademarks are the property of their respective owners. © 2019, Eurotrol B.V.

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