**DE Eurotrol GAS-ISE**

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
**EN Eurotrol GAS-ISE**

**Intended purpose**

Eurotrol GAS-ISE is an assayed quality control material intended for professional use in the verification of the precision and accuracy of combined bloodgas and electrolyte analysers: pH, pO<sub>2</sub>, pCO<sub>2</sub>, cNa<sup>+</sup>, cK<sup>+</sup>, cCl<sup>-</sup>, cLi<sup>+</sup>, cCa<sup>2+</sup>, cCa, cMg, cGlu and cLac. Eurotrol GAS-ISE is also developed to verify calibration, operating temperature and other performance characteristics.

**IVD Medical Device**

Eurotrol GAS-ISE complies with the IVD Directive 98/79/EC for Medical Devices and carries the CE mark.

**Summary**

Eurotrol GAS-ISE quality control material is available at three levels (Level 1, Level 2 and Level 3), each with known concentrations. It is intended that Eurotrol GAS-ISE should be used in the periodic verification of the precision and accuracy of different combined bloodgas and electrolyte analysers.

**Reagents**

Eurotrol GAS-ISE provides three physiological relevant levels, each ampoule holding 2.5 mL solution. Eurotrol GAS-ISE controls are prepared using salts in a physiologically buffered matrix. Tonometry with predetermined levels of oxygen and carbon dioxide balanced with nitrogen and different salt concentrations provide three distinct levels for each parameter, simulating clinically significant ranges of acid-base and electrolyte balance, respiratory function, glucose and lactate concentration. Eurotrol GAS-ISE controls contain no preservatives, viscosity adjusters or other additives that might adversely effect electrode measurements.

**Storage and stability**

Eurotrol GAS-ISE can be stored at a temperature of 15 - 30 °C. Stored unopened at this temperature it is guaranteed stable as indicated until the expiry date on the ampoules and outer box. After opening the Eurotrol GAS-ISE ampoule it is stable for 30 seconds.

**Procedures**

1. Allow the ampoule to adapt to a temperature of 25 ± 1 °C.
2. Immediately before use, shake the ampoule vigorously for at least fifteen seconds to re-equilibrate the gases with the solution. When shaken, the ampoule should be held between the thumb and index finger.
3. Swirl the ampoule gently to return the solution to the bottom of the ampoule. Allow bubbles to rise between shaking and before opening the ampoule.
4. Protect fingers with gauze, tissue or gloves.
5. **Hold ampoule with the coloured dot upside. Snap off the neck of the ampoule in the opposite direction of the coloured dot. See enclosed OPC instructions.**
6. Follow within 30 seconds the sampling procedure according to the instructions in the operating manual of your bloodgas and electrolyte analyser.
7. Please refer to local guidelines for recommended frequency of use.

**Precautions**

1. For in vitro diagnostic use only.
2. These products should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
3. Eurotrol GAS-ISE is not to be used as a calibrator.

**Expected results**

The values for each parameter at the three levels available from [www.eurotrol.com](http://www.eurotrol.com). The values are derived from randomly selected ampoules from each lot, incubated at 25 ± 1 °C and measured at 37 °C. The values are mean ranges, mean value ± 3SD. These values are a guide for evaluating the precision and accuracy of the analyser.

**Please Note**

- The values in the table are applicable only to the assigned lot number and at sea level.
- Incorrect sampling, storage, etc. may cause the readings to deviate from the target values.
- The pO<sub>2</sub> values of Eurotrol GAS-ISE vary inversely with temperature changes. For each °C change in temperatures between 18 and 26 °C, the pO<sub>2</sub> changes 1,5% in the low level, 1,1% in the normal and 0,5% in the high level. To obtain a high degree of correlation with the values in the table, the ampoules should be equilibrated as close to 25 °C as possible.
- Room temperature has no significant effect on pH values. (when opening the ampoule)
- Eurotrol GAS-ISE is sensitive to most of the instrument related factors that can cause unexpected analytical deviations.
- Eurotrol GAS-ISE controls do not contain protein and therefore may not detect malfunctions sensitive to these components.
- Eurotrol GAS-ISE is very sensitive to room air contamination. This is caused by the low solubility of oxygen in aqueous solutions which results in low buffering capacity. After 60 seconds opened, controls show 5% deviation at a low pO<sub>2</sub>, 1,1% deviation at a normal pO<sub>2</sub> and 0,4% deviation at a high pO<sub>2</sub>. Compared to oxygen, the solubility of carbon dioxide in aqueous solution is much higher and the concentration of carbon dioxide in room air is much lower. The error is less than 2% over the entire range in all three levels (20 – 70 mmHg). The effect of room air on pH is negligible.


This product has been manufactured according to Eurotrol specifications.


**For further information please contact:**

Eurotrol B.V., The Netherlands, Phone: +31 318 695 777, Fax: +31 318 695 770, E-mail: [office@eurotrol.com](mailto:office@eurotrol.com), [www.eurotrol.com](http://www.eurotrol.com)

 Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands

**Symbols used**

-  Attention, see instructions for use
-  Use by
-  In Vitro Diagnostic Medical Device
-  Manufacturer

-  Batch code
-  Temperature limitation
-  Reference number
-  CE Mark

**OPC leaflet**