DE Eurotrol Hct Control - Level B

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EN Eurotrol Hct Control - Level B

Intended Purpose

Eurotrol Hct Control - Level B is an assayed hematocrit reference material, to verify the precision and accuracy of the epoc® Blood Analysis System for the measurement of hematocrit.

IVD Medical Device

Eurotrol Hct Control - Level B complies with the IVD Directive 98/79/EC for Medical Devices and carries the CE mark. Eurotrol Hct Control - Level B complies with the following US codes of Federal Regulations (CFR): 42 CFR part 72 and 21 CFR parts 606, 640 and 820. Eurotrol Hct Control - Level B is for in vitro diagnostic use only.

Summary

Eurotrol Hct Control - Level B is an assayed hematocrit reference material, to verify the precision and accuracy of the epoc® Blood Analysis System. Each ampule containing 2.5 mL solution.

Reagents

Eurotrol Hct Control - Level B provides physiological relevant Hct levels. Eurotrol Hct Control - Level B is prepared using pure chemicals in a physiologically buffered matrix. Eurotrol Hct Control - Level B controls contain no preservatives, or other additives that might adversely effect electrode measurements.

Storage and stability

Eurotrol Hct Control - Level B can be stored at a temperature of 2-30 °C (35- 86°F). Stored unopened at this temperature product stability is guaranteed until the expiry date on the ampules and outer box. After opening the ampule, Eurotrol Hct Control - Level B is stable for 15 minutes.

Procedures

- Remove the appropriate amount of ampules from box and tray. Allow the ampule(s) to adapt to a temperature of 15-30°C (59-86 °F).
- Immediately before use, shake the ampule vigorously for at least fifteen seconds to homogenize the solution. When shaking, the ampule should be held between the thumb and index finger.
- Swirl the ampule gently to return the solution to the bottom of the ampule. Allow any bubbles to rise before opening the ampule.
- Protect fingers with gauze, tissue or gloves
- Hold ampule with non colored dot facing up. Snap off the neck of the ampule in the opposite direction of the colored dot. Please read the OPC-instruction leaflet for more details.**
- Follow the sampling procedure according to the instructions in the operating manual of your epoc Blood Analysis System.
- Please refer to local guidelines for recommended frequency of use.

	Manufacturer		CE CE mark
	Date of Manufacture		

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NL Eurotrol Hct Control - Level B

Limitations

- For in vitro diagnostic use only.
- Consult local environmental authorities for proper disposal.
- Eurotrol Hct Control - Level B is not to be used as a calibrator.

Expected Results

Values are obtained by measurements on epoc Blood Analysis System using BGEM test cards. The assigned values are available from <http://www.epocal.com>, tabulated for each fluid lot number and sensor configuration. Eurotrol Hct Control - Level B is intended for professional use in the verification of the instrument performance.

Please Note





- The values in the table are applicable only to the assigned fluid lot number, the applicable epoc Blood Analysis System sensor configuration.
- Incorect sampling, storage, or other mishandling may cause the readings to deviate from the target values.

This product has been manufactured according to Eurotrol specifications.

	Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands		
			

For further information please contact:

Eurotrol B.V. The Netherlands Phone +31 318 695 777 Fax +31 318 695 770 E-mail office@eurotrol.com www.eurotrol.com	Eurotrol Inc. USA Toll free phone 1-866-234-5754 Main phone 1-978-598-3779 Fax 1-978-598-3780 E-mail officeUSA@eurotrol.com www.eurotrol.com
--	---

	Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands		
	Attention, see instructions for use		Batch code
	Use by		Temperature limitation
	In Vitro Diagnostic Medical Device		Reference number
	Manufacturer		CE CE mark
	Date of Manufacture		

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Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, Países Baixos