Eurotrol

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### EN epoc Hematocrit Verification Fluids

**Intended Purpose** epoc® Hematocrit Verification Fluids is an assayed hematocrit reference material, to verify the precision and accuracy of the epoc® Blood Analysis System for the measurement of hematocrit.

**IVD Medical Device** epoc Hematocrit Verification Fluids complies with the IVD Directive 98/79/EC for Medical Devices and carries the CE mark. epoc Hematocrit Verification Fluids complies with the following US codes of Federal Regulations (CFR): 42 CFR part 72 and 21 CFR parts 606, 640 and 820. epoc Hematocrit Verification Fluids is for in vitro diagnostic use only.

**Summary** epoc Hematocrit Verification Fluids is a quality control material available in different levels, each with known Hct values. epoc Hematocrit Verification Fluids is an assayed hematocrit reference material, to verify the precision and accuracy of the epoc® Blood Analysis System. Each ampule containing 2.5 mL solution

**Reagents** epoc Hematocrit Verification Fluids provides physiological relevant Hct levels. epoc Hematocrit Verification Fluids is prepared using pure chemicals in a physiologically buffered matrix. Different concentrations provide five distinct Hct levels, simulating clinically significant ranges of Hematocrit. epoc Hematocrit Verification Fluids controls contain no preservatives, or other additives that might adversely effect electrode measurements.

**Storage and stability** epoc Hematocrit Verification Fluids can be stored at a temperature of 2-30 °C (35- 86°F). Stored unopened at this temperature product stability is guaranteed until the expiry date on the ampules and outer box. After opening the ampule, epoc Hematocrit Verification Fluids is stable for 15 minutes.

**Procedures**
1. Remove the appropriate amount of ampules from box and tray. Allow the ampule(s) to adapt to a temperature of 15-30 °C (59-86 °F).
2. Immediately before use, shake the ampule vigorously for at least fifteen seconds to homogenize the solution. When shaking, the ampule should be held between the thumb and index finger.
3. Swirl the ampule gently to return the solution to the bottom of the ampule. Allow any bubbles to rise before opening the ampule.
4. Protect fingers with gauze, tissue or gloves.
5. **Hold ampule with the colored dot facing up. Snap off the neck of the ampule in the opposite direction of the colored dot. Please read the OPC-instruction leaflet for more details.**
6. Follow the sampling procedure according to the instructions in the operating manual of your epoc Blood Analysis System.
7. Please refer to local guidelines for recommended frequency of use.

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**Eurotrol B.V.**