Eurotrol

Your Global Reference Point for Quality Control

HemoTrol® WB High

REF
08718734960966
CE
IVD

EN Intended Use

HemoTrol® WB is an assayed quality control material for professional use to verify the performance characteristics of the HemoCue® Hb 301 and the HemoCue® Hb 801 System. HemoTrol® WB is intended for the quantitative determination of hemoglobin.

Summary and Principle

HemoTrol® WB is a quality control material containing a known hemoglobin concentration in three clinically relevant ranges. When measured like a patient sample, the regular use of the quality control material can provide proof for proper functioning of the HemoCue® Hb 301 and the HemoCue® Hb 801 System.

Quality control materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

Reagents

Bovine blood, stabilizers, preservatives.

Precautions and Warnings

- For in vitro diagnostic use.
- Only use when the quality control has a brownish color.
- Do not freeze.
- Only mix by hand.
- Caution: HemoTrol® WB contains material of bovine origin. The product is free from transmissible spongiform encephalopathies (TSEs). However, no test method can offer complete assurance that products derived from bovine source material are free from infectious agents and may transmit infectious agents.
- Follow good laboratory practice while using the product.
- Disposition of all waste material should be in accordance with the local guidelines.
- Caution: US federal law restricts this device to sale by or on the order of a licensed practitioner.

Storage and Stability

For the shelf life at 2 to 8 °C see the expiration date on the vial label. After opening, the product is stable at 2 to 8 °C for 30 days when handled according to procedure.

Procedure

For optimal performance the product should be treated as a patient sample and instructions for use of the HemoCue® Hb 301 or HemoCue® Hb 801 System should be followed.

- Equilibrate the vial for 15 minutes to a temperature of 15 to 30 °C.
- Directly before use, shake the vial by hand for at least 30 seconds to completely mix the contents.
- Hold the vial in an upright position and tap the cap to remove liquid from the dropper tip (Figure 1.).
- Open the cap carefully without squeezing the vial (Figure 2.).
- Discard the first drop of the quality control liquid.
- Apply a drop on a hydrophobic surface and use it immediately for measurement as indicated in the analyzer user manual.
- After sampling, remove any liquid of the quality control from the colored cap with a clean tissue and recap properly.
- Return to indicated storage conditions after use. The entire procedure has to be repeated every time upon removing the quality control from storage during the entire open-vial time period.

Limitations

- The quality control material should not be used past open-vial time period or expiration date.
- The quality control material is not intended for use as calibrator / standard.
- Incomplete resuspension of the quality control material prior to use invalidates the measured sample and any following measurements made with the remaining quality control material.
- The ranges of the product are intended only as guidelines and laboratories should determine the ranges based on their own test system and tolerance limits.

Expected Values

For expected values see the batch code specific table. The means and ranges were derived from replicate analyses on the HemoCue® Hb 301 and HemoCue® Hb 801 System and are specific for that batch code of the product.

For further information or technical assistance, please contact:

NL Manufacturer Eurotrol B.V. Keplerlaan 20, 6716 BS Ede The Netherlands Phone: +31 318 695 777 E-mail: office@eurotrol.com www.eurotrol.com	Distributor Eurotrol, Inc. 850 N Black Branch Rd, Elizabethtown, KY USA Phone: +1 (502) 501-1180 E-mail: officeUSA@eurotrol.com
--	---

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state where the user and/or the patient is established.

Revision history

Revision	Effective Date	Change to previous revision
01	26-06-2018	First revision
02	08-03-2019	Addition of Rx Only symbol, the precaution that the device is intended for prescription use only and address of USA distributor.

All trademarks are the property of their respective owners. © 2018, Eurotrol B.V.

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Todas las marcas comerciales são propriedade dos respectivos proprietários. © 2018, Eurotrol B.V.