# Eurotrol

Your Global Reference Point for Quality Control

## Plasma/Low Hb Control

<span><span></span></span>	<b>Eurotrol Plasma/Low Hb Control</b>
----------------------------	---------------------------------------

**Intended purpose/Intended use**

Eurotrol Plasma/Low Hb Control is an assayed hematology control intended for professional use in the verification of the precision and accuracy of the HemoCue® Plasma/Low Hb System.

**IVD Medical Device**

The Eurotrol Plasma/Low Hb Control complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark. Eurotrol Plasma/Low Hb Control is for in vitro diagnostic use only.

**Summary**

Eurotrol Plasma/Low Hb Control is a quality control material with a known hemoglobin concentration. Eurotrol Plasma/Low Hb Control may be used in the periodic verification of the precision and accuracy of the HemoCue Plasma/Low Hb system. Please refer to local guidelines for recommended frequency of use.

**Reagents**

Eurotrol Plasma/Low Hb Control is produced at three physiologically relevant levels. Each vial containing 1.0 mL of purified bovine hemolytate.

**Storage and stability**

Eurotrol Plasma/Low Hb Control should be stored unopened in the refrigerator at 2-8 °C (36-46 °F). Stored unopened at this temperature it is guaranteed stable as indicated until the expiration date on the outer box. After opening the vial, it is stable for 30 days when properly recapped and stored at 2-8°C (36-46 °F) or at room temperature 15-30 °C (59-86 °F).

**Procedures/instructions for use**

- Allow the vial to stand for 15 minutes at room temperature 15-30 °C (59-86 °F).
- Gently mix the vial 8-10 times or until a homogeneous solution is attained, before sampling.
- Do not fill the cuvette from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the filled cuvette in the instrument for analysis within 60 seconds after filling.
- Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
- If the control does not perform as expected, review the instructions for use of the instrument to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the control still does not perform as expected, contact technical assistance.

**Precautions**

- For in vitro 1. diagnostic use only.
- Caution:** Eurotrol Plasma/Low Hb Control is an Animal Blood Product. Bovine based materials do not carry biohazards for man, such as Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). This product is free from TSE.
- It is good laboratory practice to follow the same precautions as with patient samples when handling or disposing vials.
- This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
- Eurotrol Plasma/Low Hb Control is not to be used as a calibrator.

**Expected results**

The mean assay value provided on the enclosed assay sheet is derived from replicate analyses on factory calibrated HemoCue Plasma/Low Hb systems. The values are only valid on the HemoCue Plasma/Low Hb system. Results obtained should fall within the expected range specified for the control.

For further information or technical assistance, please contact Eurotrol B.V.:

Phone +31 318 695 777
Fax +31 318 695 770
E-mail office@eurotrol.com
www.eurotrol.com

<span><span></span></span>	Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands
----------------------------	--

**ASSAYED CONTROLS FOR IN VITRO DIAGNOSTIC USE**

<span><span></span></span>	<b>Eurotrol Plasma/Low Hb Control</b>
----------------------------	---------------------------------------

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Ha bővebb tájékoztatásra vagy technikai segítségre van szüksége, forduljon az Eurotrol B.V. munkatársához: