(FR) Eurotrol Plasma/Low Hb Control

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EN Eurotrol Plasma/Low Hb Control

Intended purpose/Intended use

Eurotrol Plasma/Low Hb Control is an assayed hematology control intended for professional use in the verification of the precision and accuracy of the HemoCue® Plasma/Low Hb System.

IVD Medical Device

The Eurotrol Plasma/Low Hb Control complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark. Eurotrol Plasma/Low Hb Control complies with the following US Codes of Federal Regulations (CFR): 42 CFR part 72 and 21 CFR parts 606, 640 and 820. Eurotrol Plasma/Low Hb Control is for in vitro diagnostic use only.

Summary

Eurotrol Plasma/Low Hb Control is a quality control material with a known haemoglobin concentration. Eurotrol Plasma/Low Hb Control may be used in the periodic verification of the precision and accuracy of the HemoCue Plasma/Low Hb system. Please refer to local guidelines for recommended frequency of use.

Reagents

Eurotrol Plasma/Low Hb Control is produced at three physiologically relevant levels. Each vial containing 1.0 mL of purified bovine hemolysate.

Storage and stability

Eurotrol Plasma/Low Hb Control should be stored unopened in the refrigerator at 2-8 °C (36-46 °F). Stored unopened at this temperature it is guaranteed stable as indicated until the expiration date on the outer box. After opening the vial, it is stable for one month when properly recapped and stored at 2-30°C (36-86 °F).

Procedures/instructions for use

1. Allow the vial to stand for 15 minutes at room temperature 15-30 °C (59-86 °F).
2. Gently mix the vial 8-10 times or until a homogeneous solution is attained, before sampling.
3. Do not fill the cuvette from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the filled cuvette in the instrument for analysis within 60 seconds after filling.
4. Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
5. If the control does not perform as expected, review the instructions for use of the instrument to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the control still does not perform as expected, contact technical assistance.

Precautions

1. For in vitro diagnostic use only.
2. Caution: Eurotrol Plasma/Low Hb Control is an Animal Blood Product. Bovine based materials do not carry biohazards for man, such as Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). This product is free from TSE.
3. It is good laboratory practice to follow the same precautions as with patient samples when handling or disposing vials.
4. This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
5. Eurotrol Plasma/Low Hb Control is not to be used as a calibrator.

Expected results

The mean assay value provided on the enclosed assay sheet is derived from replicate analyses on factory calibrated HemoCue Plasma/Low Hb systems. The values are only valid on the HemoCue Plasma/Low Hb system. Results obtained should fall within the expected range specified for the control.

For further information or technical assistance, please contact:

Eurotrol B.V., The Netherlands	Eurotrol Inc., USA
Phone +31 318 695 777	Toll free phone +1 866 234 5754
	Main phone +1 978 598 3779
Fax +31 318 695 770	Fax +1 978 598 3780
E-mail office@eurotrol.com	E-mail officeUSA@eurotrol.com
www.eurotrol.com	www.eurotrol.com

 Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, the Netherlands

Symbols used

	Caution		Batch code
	Use-by date		Temperature limit
	In Vitro Diagnostic Medical Device		Catalogue number
	Manufacturer		CE Mark