

CueSee® Hyperbaric

i-STAT CLEW A37/C37/H37/J37

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	781	669 - 853	mmHg
	19 °C	66 °F	768	658 - 842	
	20 °C	68 °F	754	648 - 832	
	21 °C	70 °F	741	637 - 821	
	22 °C	72 °F	727	626 - 810	
	23 °C	73 °F	714	616 - 800	
	24 °C	75 °F	700	605 - 789	
	25 °C	77 °F	687	595 - 779	
	26 °C	79 °F	673	584 - 768	
	27 °C	81 °F	660	573 - 757	
	28 °C	82 °F	646	563 - 747	
	29 °C	84 °F	633	552 - 736	
	30 °C	86 °F	619	541 - 725	

i-STAT CLEW A37/C37/H37/J37

Cartridges types: G3+ and CG4+
Cartridge lots: D

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	752	672 - 856	mmHg
	19 °C	66 °F	743	661 - 845	
	20 °C	68 °F	734	650 - 834	
	21 °C	70 °F	725	640 - 824	
	22 °C	72 °F	716	629 - 813	
	23 °C	73 °F	707	619 - 803	
	24 °C	75 °F	698	608 - 792	
	25 °C	77 °F	689	597 - 781	
	26 °C	79 °F	680	587 - 771	
	27 °C	81 °F	672	576 - 760	
	28 °C	82 °F	663	565 - 749	
	29 °C	84 °F	654	555 - 739	
	30 °C	86 °F	645	544 - 728	

LOT 164-2-B604



Expiry date  2019/01