

## CueSee® Hyperbaric

### i-STAT CLEW A37/C37/H37/J37

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+  
Cartridge lots: A,K,L,M,N,W or Y

	Temperature		Mean	Range	Unit
pO <sub>2</sub>	18 °C	64 °F	762	665 - 849	mmHg
	19 °C	66 °F	751	655 - 839	
	20 °C	68 °F	740	644 - 828	
	21 °C	70 °F	728	634 - 818	
	22 °C	72 °F	717	623 - 807	
	23 °C	73 °F	706	612 - 796	
	24 °C	75 °F	694	602 - 786	
	25 °C	77 °F	683	591 - 775	
	26 °C	79 °F	672	580 - 764	
	27 °C	81 °F	660	570 - 754	
	28 °C	82 °F	649	559 - 743	
	29 °C	84 °F	638	549 - 733	
30 °C	86 °F	627	538 - 722		

### i-STAT CLEW A37/C37/H37/J37

Cartridges types: G3+ and CG4+  
Cartridge lots: D

	Temperature		Mean	Range	Unit
pO <sub>2</sub>	18 °C	64 °F	754	674 - 858	mmHg
	19 °C	66 °F	745	663 - 847	
	20 °C	68 °F	736	653 - 837	
	21 °C	70 °F	727	642 - 826	
	22 °C	72 °F	718	632 - 816	
	23 °C	73 °F	709	621 - 805	
	24 °C	75 °F	700	610 - 794	
	25 °C	77 °F	692	600 - 784	
	26 °C	79 °F	683	589 - 773	
	27 °C	81 °F	674	578 - 762	
	28 °C	82 °F	665	568 - 752	
	29 °C	84 °F	656	557 - 741	
30 °C	86 °F	647	547 - 731		

**LOT** 164-2-B649



Expiry date  2019/12