

## CueSee® Hyperbaric

### i-STAT CLEW A37/C37/H37/J37

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+  
Cartridge lots: A,K,L,M,N,W or Y

pO <sub>2</sub>	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	740	673 - 857	mmHg
	19 °C	66 °F	733	662 - 846	
	20 °C	68 °F	726	652 - 836	
	21 °C	70 °F	719	641 - 825	
	22 °C	72 °F	712	630 - 814	
	23 °C	73 °F	705	620 - 804	
	24 °C	75 °F	697	609 - 793	
	25 °C	77 °F	690	598 - 782	
	26 °C	79 °F	683	588 - 772	
	27 °C	81 °F	676	577 - 761	
	28 °C	82 °F	669	567 - 751	
	29 °C	84 °F	662	556 - 740	
	30 °C	86 °F	655	545 - 729	

### i-STAT CLEW A37/C37/H37/J37

Cartridges types: G3+ and CG4+  
Cartridge lots: D

pO <sub>2</sub>	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	743	681 - 865	mmHg
	19 °C	66 °F	737	671 - 855	
	20 °C	68 °F	730	660 - 844	
	21 °C	70 °F	724	649 - 833	
	22 °C	72 °F	718	639 - 823	
	23 °C	73 °F	712	628 - 812	
	24 °C	75 °F	705	618 - 802	
	25 °C	77 °F	699	607 - 791	
	26 °C	79 °F	693	596 - 780	
	27 °C	81 °F	686	586 - 770	
	28 °C	82 °F	680	575 - 759	
	29 °C	84 °F	674	564 - 748	
	30 °C	86 °F	668	554 - 738	

**LOT** 164-2-B732



Expiry date  2020/08