

CueSee® Hyperbaric

i-STAT CLEW A37/C37/H37/J37

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	753	689 - 873	mmHg
	19 °C	66 °F	746	678 - 862	
	20 °C	68 °F	740	667 - 851	
	21 °C	70 °F	733	657 - 841	
	22 °C	72 °F	726	646 - 830	
	23 °C	73 °F	720	635 - 819	
	24 °C	75 °F	713	625 - 809	
	25 °C	77 °F	706	614 - 798	
	26 °C	79 °F	700	604 - 788	
	27 °C	81 °F	693	593 - 777	
	28 °C	82 °F	686	582 - 766	
	29 °C	84 °F	680	572 - 756	
	30 °C	86 °F	673	561 - 745	

i-STAT CLEW A37/C37/H37/J37

Cartridges types: G3+ and CG4+
Cartridge lots: D

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	767	685 - 869	mmHg
	19 °C	66 °F	758	675 - 859	
	20 °C	68 °F	749	664 - 848	
	21 °C	70 °F	740	653 - 837	
	22 °C	72 °F	730	643 - 827	
	23 °C	73 °F	721	632 - 816	
	24 °C	75 °F	712	622 - 806	
	25 °C	77 °F	703	611 - 795	
	26 °C	79 °F	694	600 - 784	
	27 °C	81 °F	685	590 - 774	
	28 °C	82 °F	675	579 - 763	
	29 °C	84 °F	666	568 - 752	
	30 °C	86 °F	657	558 - 742	

LOT 164-2-B743



Expiry date  2020/10