

CueSee® Hyperbaric

i-STAT CLEW A37/C37/H37/J37

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	730	644 - 828	mmHg
	19 °C	66 °F	721	634 - 818	
	20 °C	68 °F	711	623 - 807	
	21 °C	70 °F	701	612 - 796	
	22 °C	72 °F	691	602 - 786	
	23 °C	73 °F	682	591 - 775	
	24 °C	75 °F	672	581 - 765	
	25 °C	77 °F	662	570 - 754	
	26 °C	79 °F	652	559 - 743	
	27 °C	81 °F	642	549 - 733	
	28 °C	82 °F	633	538 - 722	
	29 °C	84 °F	623	527 - 711	
	30 °C	86 °F	613	517 - 701	

i-STAT CLEW A37/C37/H37/J37

Cartridges types: G3+ and CG4+
Cartridge lots: D

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	725	658 - 842	mmHg
	19 °C	66 °F	718	648 - 832	
	20 °C	68 °F	711	637 - 821	
	21 °C	70 °F	704	626 - 810	
	22 °C	72 °F	697	616 - 800	
	23 °C	73 °F	690	605 - 789	
	24 °C	75 °F	683	595 - 779	
	25 °C	77 °F	676	584 - 768	
	26 °C	79 °F	669	573 - 757	
	27 °C	81 °F	662	563 - 747	
	28 °C	82 °F	655	552 - 736	
	29 °C	84 °F	648	541 - 725	
	30 °C	86 °F	641	531 - 715	

LOT 164-2-B822



Expiry date  2021/05