

CueSee® Hyperbaric

i-STAT CLEW A37/C37/H37/J37

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+

Cartridge lots: A, K, L, M, N, W or Y

	Temperature		Mean	Range	Unit
pO ₂	18 °C	64 °F	782	686 - 870	mmHg
	19 °C	66 °F	771	676 - 860	
	20 °C	68 °F	760	665 - 849	
	21 °C	70 °F	748	654 - 838	
	22 °C	72 °F	737	644 - 828	
	23 °C	73 °F	726	633 - 817	
	24 °C	75 °F	715	623 - 807	
	25 °C	77 °F	704	612 - 796	
	26 °C	79 °F	693	601 - 785	
	27 °C	81 °F	682	591 - 775	
	28 °C	82 °F	671	580 - 764	
	29 °C	84 °F	660	570 - 754	
	30 °C	86 °F	648	559 - 743	

LOT 164-2-B851



Expiry date  2021/12

MSP101103 Rev.1.0 20170602