

CueSee® Hyperbaric

i-STAT CLEW A37/C37/H37/J37

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	639	568 - 752	mmHg
	19 °C	66 °F	632	558 - 742	
	20 °C	68 °F	624	547 - 731	
	21 °C	70 °F	616	536 - 720	
	22 °C	72 °F	609	526 - 710	
	23 °C	73 °F	601	515 - 699	
	24 °C	75 °F	594	504 - 688	
	25 °C	77 °F	586	494 - 678	
	26 °C	79 °F	578	483 - 667	
	27 °C	81 °F	571	473 - 657	
	28 °C	82 °F	563	462 - 646	
	29 °C	84 °F	555	451 - 635	
30 °C	86 °F	548	441 - 625		

i-STAT CLEW A37/C37/H37/J37

Cartridges types: G3+ and CG4+
Cartridge lots: D

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	653	588 - 772	mmHg
	19 °C	66 °F	646	577 - 761	
	20 °C	68 °F	639	567 - 751	
	21 °C	70 °F	633	556 - 740	
	22 °C	72 °F	626	545 - 729	
	23 °C	73 °F	619	535 - 719	
	24 °C	75 °F	612	524 - 708	
	25 °C	77 °F	605	513 - 697	
	26 °C	79 °F	599	503 - 687	
	27 °C	81 °F	592	492 - 676	
	28 °C	82 °F	585	482 - 666	
	29 °C	84 °F	578	471 - 655	
30 °C	86 °F	572	460 - 644		

LOT 164-2-B905



Expiry date  2022/01

MSP101103 Rev.1.0 20170602