

CueSee® Hyperbaric

i-STAT CLEW A38/C38/H38/J38

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

| pO ₂ | Temperature | | Mean | Range | Unit |
|-----------------|-------------|-------|-----------|-----------|------|
| | °C | °F | | | |
| | 18 °C | 64 °F | 750 | 661 - 845 | mmHg |
| | 19 °C | 66 °F | 740 | 650 - 834 | |
| | 20 °C | 68 °F | 730 | 640 - 824 | |
| | 21 °C | 70 °F | 719 | 629 - 813 | |
| | 22 °C | 72 °F | 709 | 619 - 803 | |
| | 23 °C | 73 °F | 699 | 608 - 792 | |
| | 24 °C | 75 °F | 689 | 597 - 781 | |
| | 25 °C | 77 °F | 679 | 587 - 771 | |
| | 26 °C | 79 °F | 669 | 576 - 760 | |
| | 27 °C | 81 °F | 658 | 565 - 749 | |
| | 28 °C | 82 °F | 648 | 555 - 739 | |
| | 29 °C | 84 °F | 638 | 544 - 728 | |
| 30 °C | 86 °F | 628 | 534 - 718 | | |

i-STAT CLEW A38/C38/H38/J38

Cartridges types: G3+ and CG4+
Cartridge lots: D

| pO ₂ | Temperature | | Mean | Range | Unit |
|-----------------|-------------|-------|-----------|-----------|------|
| | °C | °F | | | |
| | 18 °C | 64 °F | 781 | 669 - 853 | mmHg |
| | 19 °C | 66 °F | 768 | 659 - 843 | |
| | 20 °C | 68 °F | 754 | 648 - 832 | |
| | 21 °C | 70 °F | 741 | 638 - 822 | |
| | 22 °C | 72 °F | 728 | 627 - 811 | |
| | 23 °C | 73 °F | 714 | 616 - 800 | |
| | 24 °C | 75 °F | 701 | 606 - 790 | |
| | 25 °C | 77 °F | 687 | 595 - 779 | |
| | 26 °C | 79 °F | 674 | 585 - 769 | |
| | 27 °C | 81 °F | 660 | 574 - 758 | |
| | 28 °C | 82 °F | 647 | 563 - 747 | |
| | 29 °C | 84 °F | 633 | 553 - 737 | |
| 30 °C | 86 °F | 620 | 542 - 726 | | |

LOT

164-2-B649



164-2-B649

Expiry date



2019/12

MSP101103 Rev.1.0 20170602