

CueSee® Hyperbaric

i-STAT CLEW A38/C38/H38/J38

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	748	665 - 849	mmHg
	19 °C	66 °F	739	654 - 838	
	20 °C	68 °F	729	644 - 828	
	21 °C	70 °F	720	633 - 817	
	22 °C	72 °F	711	623 - 807	
	23 °C	73 °F	701	612 - 796	
	24 °C	75 °F	692	601 - 785	
	25 °C	77 °F	683	591 - 775	
	26 °C	79 °F	673	580 - 764	
	27 °C	81 °F	664	569 - 753	
	28 °C	82 °F	655	559 - 743	
	29 °C	84 °F	645	548 - 732	
30 °C	86 °F	636	538 - 722		

i-STAT CLEW A38/C38/H38/J38

Cartridges types: G3+ and CG4+
Cartridge lots: D

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	752	668 - 852	mmHg
	19 °C	66 °F	743	657 - 841	
	20 °C	68 °F	733	647 - 831	
	21 °C	70 °F	724	636 - 820	
	22 °C	72 °F	714	625 - 809	
	23 °C	73 °F	705	615 - 799	
	24 °C	75 °F	695	604 - 788	
	25 °C	77 °F	686	594 - 778	
	26 °C	79 °F	676	583 - 767	
	27 °C	81 °F	666	572 - 756	
	28 °C	82 °F	657	562 - 746	
	29 °C	84 °F	647	551 - 735	
30 °C	86 °F	638	540 - 724		

LOT 164-2-B732



Expiry date  2020/08

MSP101103 Rev.1.0 20170602