

CueSee® Hyperbaric

i-STAT CLEW A38/C38/H38/J38

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

| pO ₂ | Temperature | | Mean | Range | Unit |
|-----------------|-------------|-------|-----------|-----------|------|
| | °C | °F | | | |
| | 18 °C | 64 °F | 755 | 673 - 857 | mmHg |
| | 19 °C | 66 °F | 746 | 663 - 847 | |
| | 20 °C | 68 °F | 737 | 652 - 836 | |
| | 21 °C | 70 °F | 728 | 642 - 826 | |
| | 22 °C | 72 °F | 719 | 631 - 815 | |
| | 23 °C | 73 °F | 709 | 620 - 804 | |
| | 24 °C | 75 °F | 700 | 610 - 794 | |
| | 25 °C | 77 °F | 691 | 599 - 783 | |
| | 26 °C | 79 °F | 682 | 589 - 773 | |
| | 27 °C | 81 °F | 673 | 578 - 762 | |
| | 28 °C | 82 °F | 664 | 567 - 751 | |
| | 29 °C | 84 °F | 655 | 557 - 741 | |
| 30 °C | 86 °F | 645 | 546 - 730 | | |

i-STAT CLEW A38/C38/H38/J38

Cartridges types: G3+ and CG4+
Cartridge lots: D

| pO ₂ | Temperature | | Mean | Range | Unit |
|-----------------|-------------|-------|-----------|-----------|------|
| | °C | °F | | | |
| | 18 °C | 64 °F | 748 | 679 - 863 | mmHg |
| | 19 °C | 66 °F | 740 | 668 - 852 | |
| | 20 °C | 68 °F | 733 | 657 - 841 | |
| | 21 °C | 70 °F | 726 | 647 - 831 | |
| | 22 °C | 72 °F | 718 | 636 - 820 | |
| | 23 °C | 73 °F | 711 | 626 - 810 | |
| | 24 °C | 75 °F | 704 | 615 - 799 | |
| | 25 °C | 77 °F | 696 | 604 - 788 | |
| | 26 °C | 79 °F | 689 | 594 - 778 | |
| | 27 °C | 81 °F | 682 | 583 - 767 | |
| | 28 °C | 82 °F | 674 | 573 - 757 | |
| | 29 °C | 84 °F | 667 | 562 - 746 | |
| 30 °C | 86 °F | 660 | 551 - 735 | | |

LOT 164-2-B743



Expiry date  2020/10

MSP101103 Rev.1.0 20170602