

CueSee® Hyperbaric

i-STAT CLEW A38/C38/H38/J38

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	720	639 - 823	mmHg
	19 °C	66 °F	711	629 - 813	
	20 °C	68 °F	702	618 - 802	
	21 °C	70 °F	693	608 - 792	
	22 °C	72 °F	684	597 - 781	
	23 °C	73 °F	675	586 - 770	
	24 °C	75 °F	666	576 - 760	
	25 °C	77 °F	657	565 - 749	
	26 °C	79 °F	648	554 - 738	
	27 °C	81 °F	639	544 - 728	
	28 °C	82 °F	630	533 - 717	
	29 °C	84 °F	621	523 - 707	
	30 °C	86 °F	612	512 - 696	

i-STAT CLEW A38/C38/H38/J38

Cartridges types: G3+ and CG4+
Cartridge lots: D

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	731	650 - 834	mmHg
	19 °C	66 °F	722	640 - 824	
	20 °C	68 °F	713	629 - 813	
	21 °C	70 °F	704	618 - 802	
	22 °C	72 °F	695	608 - 792	
	23 °C	73 °F	686	597 - 781	
	24 °C	75 °F	677	587 - 771	
	25 °C	77 °F	668	576 - 760	
	26 °C	79 °F	659	565 - 749	
	27 °C	81 °F	650	555 - 739	
	28 °C	82 °F	641	544 - 728	
	29 °C	84 °F	632	533 - 717	
	30 °C	86 °F	623	523 - 707	

LOT 164-2-B822



Expiry date  2021/05

MSP101103 Rev.1.0 20170602