

CueSee® Hyperbaric

i-STAT CLEW A38/C38/H38/J38

Cartridges types: CG8+, EG7+, EG6+, G3+ and CG4+
Cartridge lots: A,K,L,M,N,W or Y

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	632	567 - 751	mmHg
	19 °C	66 °F	625	556 - 740	
	20 °C	68 °F	619	546 - 730	
	21 °C	70 °F	612	535 - 719	
	22 °C	72 °F	605	525 - 709	
	23 °C	73 °F	598	514 - 698	
	24 °C	75 °F	591	503 - 687	
	25 °C	77 °F	585	493 - 677	
	26 °C	79 °F	578	482 - 666	
	27 °C	81 °F	571	471 - 655	
	28 °C	82 °F	564	461 - 645	
	29 °C	84 °F	558	450 - 634	
30 °C	86 °F	551	440 - 624		

i-STAT CLEW A38/C38/H38/J38

Cartridges types: G3+ and CG4+
Cartridge lots: D

pO ₂	Temperature		Mean	Range	Unit
	°C	°F			
	18 °C	64 °F	655	581 - 765	mmHg
	19 °C	66 °F	647	571 - 755	
	20 °C	68 °F	639	560 - 744	
	21 °C	70 °F	631	549 - 733	
	22 °C	72 °F	623	539 - 723	
	23 °C	73 °F	615	528 - 712	
	24 °C	75 °F	607	518 - 702	
	25 °C	77 °F	599	507 - 691	
	26 °C	79 °F	591	496 - 680	
	27 °C	81 °F	583	486 - 670	
	28 °C	82 °F	575	475 - 659	
	29 °C	84 °F	567	464 - 648	
30 °C	86 °F	559	454 - 638		

LOT 164-2-B905



Expiry date  2022/01

MSP101103 Rev.1.0 20170602