

CueSee® Hypoxic

X-Low

LOT 16301812



REF 166.001.005

Radiometer ABL 835 FLEX

Analyte	Mean	Range	Unit
pO_2	20.6	8.6 – 32.6	mmHg

Abbott i-STAT® system CLEW: A37

Cartridge types: CG8+, EG7+, EG6+, G3+ & CG4+ Lots: A, K, L, M, N, W, Y

Analyte	Mean	Range	Unit
pO_2	16	5 – 28	mmHg

Cartridge types: G3+ & CG4+ Lot: D

Analyte	Mean	Range	Unit
pO_2	16	5 – 28	mmHg

epoc® Blood Analysis System

See instrument manufacturer's website for Assigned values

Assigned values are established on the instruments listed. For other instruments mean values and ranges can be established by the end user.

AN01183A01 – 20180615