

CueSee® Hypoxic

X-Low

LOT 16301911



REF 166.001.005

Radiometer ABL 835 FLEX

Analyte	Mean	Range	Unit
pO_2	18.6	6.6 – 30.6	mmHg

Abbott i-STAT® system CLEW: A37

Cartridge types: CG8+, EG7+, EG6+, G3+ & CG4+ Lots: A, K, L, M, N, W, Y

Analyte	Mean	Range	Unit
pO_2	14	5 – 26	mmHg

Cartridge types: G3+ & CG4+ Lot: D

Analyte	Mean	Range	Unit
pO_2	15	5 – 27	mmHg

epoc® Blood Analysis System

See instrument manufacturer's website for Assigned values

Assigned values are established on the instruments listed.
For other instruments mean values and ranges can be established by the end user.

AN01183A01 – 20180615