


Barcodes for Target Ranges


GlucoTrol-NG Level 3

**HemoCue B-Glucose systems
HemoCue Glucose 201 systems**

whole blood "B-Glu"

LOT 82370  2019-10


	Mean	Range
mmol/L	9.9	8.7 - 11.1
mg/dL	179	157 - 201




2 1 9 3 5 1 8 2 3 7 0 1 7 9 1 5 7 2 0 1

**HemoCue B-Glucose systems
HemoCue Glucose 201 systems**

plasma equivalent "P-Equiv"

LOT 82370  2019-10


	Mean	Range
mmol/L	11.0	9.8 - 12.2
mg/dL	199	177 - 221




2 1 9 3 5 1 8 2 3 7 0 1 9 9 1 7 7 2 2 1

HemoCue Glucose 201RT systems

whole blood "B-Glu"

LOT 82370  2019-10


	Mean	Range
mmol/L	9.6	8.4 - 10.8
mg/dL	173	151 - 195




7 1 9 3 5 1 8 2 3 7 0 2 0 0 1 5 1 1 9 5

HemoCue Glucose 201RT systems

plasma equivalent "P-Equiv"

LOT 82370  2019-10

	Mean	Range
mmol/L	10.6	9.4 - 11.8
mg/dL	192	170 - 214



7 1 9 3 5 1 8 2 3 7 0 2 0 1 1 7 0 2 1 4

AN01068A01 20180130