


## Barcodes for Target Ranges


### GlucoTrol-NG Level 3

**HemoCue B-Glucose systems  
HemoCue Glucose 201 systems**

**whole blood "B-Glu"**

LOT 73362     2019-01


	Target	Expected range
mmol/L	10.0	8.8 - 11.2
mg/dL	181	159 - 203




2 1906373362181159203

**HemoCue B-Glucose systems  
HemoCue Glucose 201 systems**

**plasma equivalent "P-Equiv"**

LOT 73362     2019-01


	Target	Expected range
mmol/L	11.1	9.9 - 12.3
mg/dL	201	179 - 223




2 1906373362201179223

**HemoCue B-Glucose systems  
HemoCue Glucose 201 systems**

**whole blood "B-Glu"**

LOT 73362     2019-01


mmol/L	$\bar{x}$	$\pm$
mg/dL	10.0	1.2
	181	22




2 1906373362181159203

**HemoCue B-Glucose systems  
HemoCue Glucose 201 systems**

**plasma equivalent "P-Equiv"**

LOT 73362     2019-01


mmol/L	$\bar{x}$	$\pm$
mg/dL	11.1	1.2
	201	22




2 1906373362201179223

**HemoCue Glucose 201RT systems**

**whole blood "B-Glu"**

LOT 73362     2019-01


	Target	Expected range
mmol/L	9.8	8.6 - 11.0
mg/dL	176	154 - 198




7 1906373362200154198

**HemoCue Glucose 201RT systems**

**plasma equivalent "P-Equiv"**

LOT 73362     2019-01


	Target	Expected range
mmol/L	10.8	9.6 - 12.0
mg/dL	195	173 - 217




7 1906373362201173217

**HemoCue Glucose 201RT systems**

**whole blood "B-Glu"**

LOT 73362     2019-01


mmol/L	$\bar{x}$	$\pm$
mg/dL	9.8	1.2
	176	22




7 1906373362200154198

**HemoCue Glucose 201RT systems**

**plasma equivalent "P-Equiv"**

LOT 73362     2019-01

mmol/L	$\bar{x}$	$\pm$
mg/dL	10.8	1.2
	195	22



7 1906373362201173217