



Barcodes for Target Ranges
 GlucoTrol-NG Level 4

HemoCue B-Glucose systems
HemoCue Glucose 201 systems

whole blood "B-Glu"

LOT 73363  2019-01


	Target	Expected range
mmol/L	17.2	14.9 - 19.5
mg/dL	310	268 - 352




3 1906373363310268352

HemoCue B-Glucose systems
HemoCue Glucose 201 systems

plasma equivalent "P-Equiv"

LOT 73363  2019-01


	Target	Expected range
mmol/L	19.1	16.8 - 21.4
mg/dL	344	302 - 386




3 1906373363344302386

HemoCue B-Glucose systems
HemoCue Glucose 201 systems

whole blood "B-Glu"

LOT 73363  2019-01


	\bar{x}	\pm
mmol/L	17.2	2.3
mg/dL	310	42




3 1906373363310268352

HemoCue B-Glucose systems
HemoCue Glucose 201 systems

plasma equivalent "P-Equiv"

LOT 73363  2019-01


	\bar{x}	\pm
mmol/L	19.1	2.3
mg/dL	344	42




3 1906373363344302386

HemoCue Glucose 201RT systems

whole blood "B-Glu"

LOT 73363  2019-01


	Target	Expected range
mmol/L	16.9	14.6 - 19.2
mg/dL	304	262 - 346




7 1906373363300262346

HemoCue Glucose 201RT systems

plasma equivalent "P-Equiv"

LOT 73363  2019-01


	Target	Expected range
mmol/L	18.7	16.4 - 21.0
mg/dL	337	295 - 379




7 1906373363301295379

HemoCue Glucose 201RT systems

whole blood "B-Glu"

LOT 73363  2019-01


	\bar{x}	\pm
mmol/L	16.9	2.3
mg/dL	304	42




7 1906373363300262346

HemoCue Glucose 201RT systems

plasma equivalent "P-Equiv"

LOT 73363  2019-01

	\bar{x}	\pm
mmol/L	18.7	2.3
mg/dL	337	42



7 1906373363301295379