


## Barcodes for Target Ranges


### Glucotrol-NG Level 2

**HemoCue B-Glucose systems  
HemoCue Glucose 201 systems**

**whole blood "B-Glu"**

**LOT** 74569  2019-04


	Target	Expected range
mmol/L	6.0	5.1 - 6.9
mg/dL	109	92 - 126




1 19 15874569109092126

**HemoCue B-Glucose systems  
HemoCue Glucose 201 systems**

**plasma equivalent "P-Equiv"**

**LOT** 74569  2019-04


	Target	Expected range
mmol/L	6.7	5.8 - 7.6
mg/dL	121	104 - 138




1 19 15874569121104138

**HemoCue Glucose 201RT systems**

**whole blood "B-Glu"**

**LOT** 74569  2019-04


	Target	Expected range
mmol/L	5.7	4.8 - 6.6
mg/dL	103	86 - 120




7 19 15874569100086120

**HemoCue Glucose 201RT systems**

**plasma equivalent "P-Equiv"**

**LOT** 74569  2019-04

	Target	Expected range
mmol/L	6.4	5.5 - 7.3
mg/dL	115	98 - 132



7 19 15874569101098132