


Barcodes for Target Ranges


Glucotrol-NG Level 1

**HemoCue B-Glucose systems
HemoCue Glucose 201 systems**

whole blood "B-Glu"

LOT 73359  2019-01


	Target	Expected range
mmol/L	2.0	1.2 - 2.8
mg/dL	36	22 - 50




0 1906373359036022050

**HemoCue B-Glucose systems
HemoCue Glucose 201 systems**

plasma equivalent "P-Equiv"

LOT 73359  2019-01


	Target	Expected range
mmol/L	2.2	1.4 - 3.0
mg/dL	40	26 - 54




0 1906373359040026054

HemoCue Glucose 201RT systems

whole blood "B-Glu"

LOT 73359  2019-01


	Target	Expected range
mmol/L	2.4	1.6 - 3.2
mg/dL	43	29 - 57




7 1906373359000029057

HemoCue Glucose 201RT systems

plasma equivalent "P-Equiv"

LOT 73359  2019-01

	Target	Expected range
mmol/L	2.6	1.8 - 3.4
mg/dL	47	33 - 61



7 1906373359001033061