


Barcodes for Target Ranges


HemoTrol Low (Level 1)

HemoCue Hemoglobin systems

Hb-201+ analyzers

LOT 73301  2019-04


	Target	Expected range
g/dL	8.0	7.6 - 8.4
g/L	80	76 - 84
mmol/L	5.0	4.8 - 5.2




0 19 15 87 33 01 58 03 76 08 4

HemoCue Hemoglobin systems

Hb-201+ analyzers

LOT 73301  2019-04

	Target	±
g/dL	8.0	0.4
g/L	80	4
mmol/L	5.0	0.2



0 19 15 87 33 01 58 03 76 08 4