


Barcodes for Target Ranges


HemoTrol Normal (Level 2)

HemoCue Hemoglobin systems

Hb-201+ analyzers

	LOT	73602		2019-04
--	------------	--------------	---	----------------


	Target	Expected range
g/dL	11.9	11.3 -12.5
g/L	119	113 - 125
mmol/L	7.4	7.0 - 7.8




1 19 15 87 36 02 61 94 13 12 5

HemoCue Hemoglobin systems

Hb-201+ analyzers

	LOT	73602		2019-04
--	------------	--------------	---	----------------

	Target	±
g/dL	11.9	0.6
g/L	119	6
mmol/L	7.4	0.4



1 19 15 87 36 02 61 94 13 12 5